

Aflatoxin HPLC Column from LCTech, Leaders in Mycotoxin Analysis



The New Aflatoxin HPLC column from LCTech allows faster analysis, improving efficiency and saving solvents.

The integral guard has zero impact on resolution and helps protect the column from premature failure due to matrix contamination.


The standard method uses less than 15% Acetonitrile, but still allows rapid analysis.

Combined with the UVE Postcolumn photoderivatizer, it is the perfect complement to LCTech's AflaClean™ Immunoaffinity columns

Highlights:

- Separation in Under 8 minutes
- Reduces solvent consumption by 50%
- Saves time and improves productivity
- Low AcN in mobile phase requires <1.35ml/sample run
- Integral Guard system

Part No	Ordering Information	Price £
10522	HPLC-Column for Aflatoxin Analysis, Reversed Phase, C18, 150 x 4.6 mm	£410.00
10750	Guard Cartridge Holder	£110.00
10523	Guard Cartridges for a.m. HPLC column (3 pcs), C18, 8 x 4 mm	£180.00



www.arcsciences.com

ARC House
PO Box 275
Alton
Hampshire
GU34 9FJ

Phone > +44 (0) 1420 549922
Fax > +44 (0) 1420 84254

VAT No: 935779078
Company Reg: 6547940

enquiries@arcsciences.com